# Food Science of Animal Resources

Food Sci. Anim. Resour. 2021 January 41(1):85~94 DOI https://doi.org/10.5851/kosfa.2020.e80





# Comparison of Seven Commercial TaqMan Master Mixes and Two Real-Time PCR Platforms Regarding the Rapid Detection of Porcine DNA

Soo Ji Kang, Chan Song Jang, Ji Min Son, and Kwang Won Hong\*

Department of Food Science and Biotechnology, College of Life Science and Biotechnology, Dongguk University, Goyang 10326, Korea



June 8, 2020 Received Revised August 19, 2020 Accepted September 11, 2020

\*Corresponding author: Kwang Won Hong Department of Food Science and Biotechnology, College of Life Science and Biotechnology, Dongguk University, Goyang 10326, Korea Tel: +82-31-961-5140 Fax: +82-31-961-5108

E-mail: hkwon@dongguk.edu

#### \*ORCID

Soo Ji Kang https://orcid.org/0000-0002-1783-2116 Chan Song Jang https://orcid.org/0000-0002-8590-605X Ji Min Son https://orcid.org/0000-0001-9850-719X Kwang Won Hong  $https://orcid.org/\overline{0000-0001-8717-8036}$ 

**Abstract** A pig-specific real-time PCR assay based on the mitochondrial ND5 gene was developed to detect porcine material in food and other products. To optimize the performance of assay, seven commercial TaqMan master mixes and two real-time PCR platforms (Applied Biosystems StepOnePlus and Bio-rad CFX Connect) were used to evaluate the limit of detection (LOD) as well as the PCR efficiency and specificity. The LODs and PCR efficiencies for the seven master mixes on two platforms were 0.5-5 pg/reaction and 84.96%–108.80%, respectively. Additionally, non-specific amplifications of DNA from other animal samples (human, dog, cow, and chicken) were observed for four master mixes. These results imply that the sensitivity and specificity of a real-time PCR assay may vary depending on master mix and platform used. The best combination of master mix and real-time PCR platform can accurately detect 0.5 pg porcine DNA, with a PCR efficiency of 100.49%.

**Keywords** master mix, real-time PCR, species identification, porcine DNA

#### Introduction

Correctly identifying meat species in food products is very important for authenticating food, promoting food safety, and preventing food adulteration. Meat species identification is a critical issue because of the different forms of meat adulteration, including the replacement of expensive meat with cheaper meat, the presence of less meat than indicated on the product label, and the inclusion of meat in non-meat (vegetarian) products (Zia et al., 2020). Furthermore, accurate meat species identification is important for satisfying religious requirements for certain foods (e.g., Halal meat). Islamic law strictly forbids the consumption of some meat products, especially pork. Thus, there are authenticity problem and religious reason for protecting consumers by detecting pork in food products.

Among the various analytical methods available for detecting meat species in foods,

highly sensitive and specific DNA-based methods have been commonly applied (El Sheikha et al., 2017). PCR-based methods involving random amplified polymorphic DNA (Arslan et al., 2005; Mane et al., 2008), restriction fragment length polymorphisms (Hossain et al., 2016; Rahmati et al., 2016), DNA barcoding (Kane and Hellberg 2016; Naaum et al., 2018), and real-time PCR (Amaral et al., 2017; Mohamad et al., 2018) are frequently used because they enable rapid and precise detection of meat species. Among these options, real-time PCR with a species-specific primer and a TaqMan probe is the most suitable and widely used method for identifying meat species (Ali et al., 2012; Kim et al., 2016). Furthermore, there are numerous available commercial master mixes and platforms, with increasing interest in the TaqMan real-time PCR assay.

Several previous studies investigated the effects of different master mixes and platforms on real-time PCR performance characteristics. For example, Sohni et al. (2008) evaluated five commercial real-time PCR reagents used for detecting *Bacillus anthracis* by comparing their PCR efficiencies and limits of detection (LODs). Stephens et al. (2010) compared five master mixes used for detecting the Ebola virus regarding their performance characteristics such as sensitivity and PCR efficiency. Meanwhile, Buzard et al. (2012) conducted a multi-platform comparison of nine commercial master mixes used for detecting bioterrorism agents. Furthermore, Eischeid and Kasko (2015) compared the utility of four master mixes for identifying a shrimp allergen in a real-time PCR assay. All of these studies emphasized the importance of choosing a suitable master mix and platform, both of which can influence the sensitivity and efficiency of a PCR assay. However, there have been relatively few comparative studies regarding real-time PCR master mixes and platforms used for identifying meats, especially pork. Moreover, most of these studies did not test whether the PCR assay was specific for meat from a particular animal species.

Therefore, the main objective of this study was to optimize a real-time PCR assay for detecting porcine material. Seven commercial master mixes were compared regarding specific performance criteria, including LOD, PCR efficiency, specificity, total cost, and time. The reliability of the data was confirmed with two real-time PCR platforms.

#### **Materials and Methods**

## **Commercial master mixes**

The following seven commercial TaqMan master mixes were evaluated: TaqMan Universal PCR Master Mix (Applied Biosystems, Waltham, MA, USA), MG 2X qPCR MasterMix (TaqMan) with ROX (CancerROP, Seoul, Korea), Express qPCR Supermix Universal (Invitrogen, Waltham, MA, USA), PowerAmp Real-time PCR Master Mix II (Kogene Biotech, Seoul, Korea), Luna Universal Probe qPCR Master Mix (New England Biolabs, Ipswich, MA, USA), QuantiNova Probe PCR Kit (Qiagen, Hilden, Germany), and Premix Ex Taq (Probe qPCR), ROX plus (Takara, Shiga, Japan).

#### Sample collections

Four raw meat samples including pig (Sus scrofa domesticus), dog (Canis familiaris), chicken (Gallus gallus), and cow (Bos taurus) were purchased from local supermarkets in Korea, and human cheek cells were obtained after rinsing the mouth with 1 mL 8% NaCl.

Meat-processed foods including 4 types of pork-containing products (dumpling, ham, pork cutlet, and sausage), 3 types of beef-containing products (beef curry, beef stock, and meatballs), and 3 types of chicken-containing products (chicken teriyaki, chicken sausage, and chicken stock) were purchased from local markets in Korea.

#### **DNA** extraction

Genomic DNA was extracted from 200 mg of finely ground samples. For all samples, DNA was extracted with the PowerPrep<sup>TM</sup> DNA Extraction from Food and Feed Kit (Kogene Biotech) according to the manufacturer's instruction. The concentration and quality of extracted DNA were determined by Nanodrop 2000 spectrophotometer (Thermo Fisher Scientific, Waltham, MA, USA). The DNA concentration of all samples was adjusted to 50 ng/μL and a 10-fold serial dilution series of porcine DNA was prepared to generate a real-time PCR standard curve.

#### Primer and probe design

The porcine-specific primer set [5'-CGCCTCACTCACATTAACCA-3' (forward) and 5'- AAGGGGACTAGGCTGAGA GTG-3' (reverse)] and TaqMan probe [5'-FAM- CACTGACTATTCTAACCATCCCAA-BHQ1-3'] were designed as follows. Various porcine DNA sequences from the National Center for Biotechnology Information GenBank database (http://www.ncbi.nlm.nih.gov/Genbank) were aligned with the ClustalW program (http://www.ebi.ac.uk/clustalw/). The primer set was designed to specifically target the conserved region of the *ND5* gene (NADH dehydrogenase subunit 5) from the porcine mitochondrial genome to produce a 141-bp amplicon. The TaqMan probe was tag with FAM and BHQ1 at the 5' and 3' ends, respectively. The primer set and probe were synthesized by Bioneer (Deajeon, Korea).

#### Real-time PCR assay and data analysis

The reaction mixture for the real-time PCR assay comprised 10 μL each master mix, 500 nM primer set, 500 nM TaqMan probe, 1 μL DNA (10-fold serial dilution series), and distilled water for a final volume of 20 μL. The manufacturer's recommended thermal cycling conditions used in this study are listed in Table 1. To determine the LODs and PCR efficiencies, all samples of the 10-fold serial dilution series were analyzed in triplicate on the same 96-well optical reaction plate (Applied Biosystems). The real-time PCR assay was completed using the StepOnePlus<sup>TM</sup> Real-Time PCR system (Applied Biosystems) and the CFX Connect<sup>TM</sup> Real-Time PCR System (Bio-Rad, Hercules, CA, USA).

All data were analyzed with StepOnePlus<sup>TM</sup> Software (version 2.3) (Applied Biosystems) and CFX Manager<sup>TM</sup> Software (Bio-Rad). Moreover, these programs automatically generated the standard curve and determined the PCR efficiency, which was calculated as  $E=-1+10^{(-1/\text{slope})}$ .

Table 1. Summary of the seven commercial master mixes evaluated in this study

Manufacturer	Master mix	Manufacturer's recommended thermal cycling conditions
Applied Biosystems	TaqMan Universal PCR Master Mix	50°C for 2 min, 95°C for 10 min, and 35 cycles of 95°C for 15 s and 60°C for 1 min
CancerROP	MG 2X qPCR MasterMix (TaqMan) with ROX	95°C for 5 min and 35 cycles of 95°C for 30 s and 60°C for 1 min
Invitrogen	Express qPCR Supermix Universal	$50^{\circ}\text{C}$ for 2 min, 95°C for 2 min, and 35 cycles of 95°C for 15 s and 60°C for 1 min
Kogene Biotech	PowerAmp Real-time PCR Master Mix II	$50^{\circ}\text{C}$ for 2 min, 95°C for 10 min, and 35 cycles of 95°C for 15 s and 60°C for 1 min
New England Biolabs	Luna Universal Probe qPCR Master Mix	$95^{\circ}\text{C}$ for $60\text{ s}$ and $35\text{ cycles}$ of $95^{\circ}\text{C}$ for $15\text{ s}$ and $60^{\circ}\text{C}$ for $30\text{ s}$
Qiagen	QuantiNova Probe PCR Kit	$95^{\circ}\text{C}$ for 2 min and 35 cycles of $95^{\circ}\text{C}$ for 5 s and $60^{\circ}\text{C}$ for $10~\text{s}$
Takara	Premix Ex Taq (Probe qPCR), ROX plus	$95^{\circ}\text{C}$ for 20 s and 35 cycles of $95^{\circ}\text{C}$ for 1 s and $60^{\circ}\text{C}$ for 20 s

#### **Results and Discussion**

## Limit of detection (LOD) and PCR efficiency

To compare the LODs and PCR efficiencies, 10-fold serial dilutions of porcine DNA (0.0005–50 ng/μL) were analyzed in a real-time PCR assay involving seven master mixes and two platforms (Table 2). All master mixes performed stably on both platforms, but the LODs of the master mixes ranged from 0.5 to 5 pg/reaction (rxn). Sensitive and reliable amplifications were observed for the Kogene Biotech, Invitrogen, Qiagen, and New England Biolabs master mixes, with an LOD of 0.5 pg/rxn on both platforms. In contrast, the LOD for the Applied Biosystems and CancerROP master mixes differed depending on the platform. Specifically, the LOD for the Applied Biosystems master mix was 10 times lower on the StepOnePlus platform (0.5 pg/rxn) than on the CFX Connect platform (5 pg/rxn). Conversely, the LOD of the CancerROP master mix was 10 times lower on the CFX Connect platform (0.5 pg/rxn) than on the StepOnePlus platform (5 pg/rxn). The Takara master mix was the least sensitive, with an LOD of 5 pg/rxn. The PCR efficiencies of the seven master mixes ranged from 84.96% to 108.80% depending on the master mixes and platforms (Fig. 1). The Kogene Biotech master mix on the CFX Connect platform performed best, with an efficiency of 100.49% (correlation coefficient, r²=0.9997). Meanwhile, the worst performance was observed for the Applied Biosystems master mix on the CFX Connect platform, with an efficiency of 84.96% (correlation coefficient, r²=0.9943).

The data presented herein revealed that the sensitivity and efficiency of a real-time PCR assay varied depending on the master mix and platform used. The observed differences were due to the DNA polymerase and buffer in the master mix, both of which influence the amplification efficiency and ability to detect specific DNA sequences (Wolffs et al., 2004). Previous studies that evaluated the utility of DNA polymerases for amplifying DNA samples collected during forensic analyses

Table 2. Comparison of the sensitivity of seven commercial master mixes for detecting porcine DNA by real-time PCR

Manufacturer	Platform	C <sub>T</sub> values for different concentrations of porcine DNA (ng/μL)					
Manufacturer	Piatiorm -	5×10 <sup>1</sup>	5×10 <sup>0</sup>	5×10 <sup>-1</sup>	5×10 <sup>-2</sup>	5×10 <sup>-3</sup>	5×10 <sup>-4</sup>
Applied Biosystems	StepOnePlus	17.77±0.53	$20.60\pm0.47$	$23.97 \pm 0.46$	$28.44 \pm 0.93$	$31.00 \pm 0.08$	34.79±0.33
_	CFX Connect	$17.38 \pm 0.23$	$20.01 \pm 0.18$	$24.65 \pm 0.04$	$27.92 \pm 0.32$	$32.15 \pm 0.04$	ND
CancerROP	StepOnePlus	18.54±0.09	20.56±0.12	25.44±0.47	29.18±0.48	32.09±0.72	ND
	CFX Connect	$15.72 \pm 0.44$	$18.23 \pm 0.64$	$22.54 \pm 0.52$	$26.52 \pm 0.02$	$29.91 \pm 0.17$	$33.41 \pm 0.52$
Invitrogen	StepOnePlus	17.01±0.16	19.79±0.34	23.71±0.21	26.91±0.22	30.86±0.30	33.78±0.21
	CFX Connect	$17.60 \pm 0.34$	$21.05 \pm 0.02$	$24.57 \pm 0.13$	$28.45 \pm 0.13$	$31.40 \pm 0.04$	$34.75 \pm 0.13$
Kogene Biotech	StepOnePlus	15.35±0.37	18.59±0.35	21.96±0.64	25.13±0.15	28.92±0.05	32.01±0.08
_	CFX Connect	$17.40 \pm 0.45$	20.55±0.42	23.97±0.41	$27.20 \pm 0.03$	$30.75 \pm 0.12$	33.81±0.43
New England Biolabs	StepOnePlus	$17.27 \pm 0.04$	$19.82 \pm 0.24$	$24.39 \pm 0.09$	$27.85 \pm 0.05$	31.16±0.19	34.38±1.03
_	CFX Connect	$17.61\pm0.11$	19.11±0.69	$22.38 \pm 0.34$	$25.89 \pm 0.48$	29.65±0.54	32.56±0.16
Qiagen	StepOnePlus	17.19±0.23	20.05±0.27	24.09±0.57	27.23±0.93	31.02±0.69	33.67±0.35
	CFX Connect	$17.08 \pm 0.26$	$20.81 \pm 0.06$	$22.81 \pm 0.13$	$26.62 \pm 0.32$	$30.36 \pm 0.45$	$33.36 \pm 0.17$
Takara	StepOnePlus	18.69±0.37	22.00±0.45	26.47±0.18	29.91±0.69	33.28±0.21	ND
	CFX Connect	16.68±0.18	19.65±0.25	22.77±0.34	26.24±0.23	29.02±0.21	ND

Average  $C_T$  values (mean $\pm$ SD) for 50, 5, 0.5, 0.05, 0.005, and 0.0005 ng/ $\mu$ L porcine DNA analyzed in triplicate. ND, not detected.




Fig. 1. Standard curves for seven master mixes (a–g) tested on the StepOnePlus and CFX Connect platforms. The real-time PCR assay was completed in triplicate using 10-fold serial dilutions of porcine DNA. Error bars are not shown because the symbol is larger than the error bar.

compared nine DNA polymerases (Ampli*Taq* Gold, Bio-X-Act Short, Ex*Taq* Hot Start, KAPA2G Robust, Omni*Taq*, PicoMaxx High Fidelity, r*Tth*, *Taq*, and *Tth*). The resulting LODs and PCR efficiencies differed depending on the DNA polymerase used (Hedman et al., 2009; Hedman et al., 2010). Other studies demonstrated that the master mix buffer components, such as Mg<sup>2+</sup> and bovine serum albumin, influence DNA polymerase activity (Bustin, 2004; Kreader, 1996). Therefore, the variability in the LODs and PCR efficiencies in the current study was because the seven tested master mixes comprise a different DNA polymerase–buffer system. Moreover, the suitability of a particular DNA polymerase–buffer system may depend on the DNA target.

This comparative study of seven master mixes may be useful for the development of a very sensitive real-time PCR assay for detecting porcine material. A comparison with several previous real-time PCR-based studies of porcine detection indicated that the LOD (0.0005 ng/rxn) of this study is lower than that (0.01 ng/rxn) of a previous study by Rodríguez et al. (2005), but is similar to the LOD (0.0001 ng/rxn) of a study by Kesmen et al. (2009), which involved the most sensitive real-time PCR system for detecting porcine material. Additionally, appropriate PCR efficiencies were reportedly between 90% and 110% (Adams, 2006). In this study, the PCR efficiencies of most of the master mix and platform combinations were between 84.96% and 108.80%. Moreover, the highest efficiency (100.49%) observed for the Kogene Biotech master mix on the CFX Connect platform was greater than the previously reported PCR efficiencies of 75.83% (Rodríguez et al., 2005), 91.57% (Sakai et al., 2011), and 103.98% (Kesmen et al., 2009).

## **Specificity test**

The specificity of the seven master mixes was tested with 50 ng/μL DNA extracted from four animal species (human, dog, cow, and chicken), with porcine DNA used as a positive control. All samples were analyzed in triplicate along with a notemplate control, and the results are presented in Table 3. The positive control was amplified in all seven master mixes, with C<sub>T</sub> values between 15.35±0.37 and 18.69±0.37. A lack of non-specific amplification was observed for the Applied Biosystems, Takara, and Kogene Biotech master mixes on both platforms. Meanwhile, non-specific amplifications were detected for the Invitrogen and New England Biolabs master mixes, but only on the CFX Connect platform, while they were observed for the Qiagen and CancerROP master mixes on both platforms.

An additional sequence alignment analysis of five animal species (pig, human, dog, cow, and chicken) with the ClustalW program was used to verify the specificity of the primer sets and probes. The primer and probe sequences were completely complementary to the mitochondrial DNA of *Sus scrofa domesticus* (NC012096.1, AF486858.1, AF486866.1, AY574046.1, DQ518915.2, EU117375.1, KJ746666.1, and KC469587.1), but not to the sequences from other species, specifically *Homo sapiens* (GU170821.1), *Canis familiaris* (AY729880.1), *Bos taurus* (GU947021.1), and *Gallus gallus* (KM096864.1). A conventional PCR was applied to confirm the specificity of the primer set developed in this study. The electrophoretic separation of amplicon revealed a lack of non-specific products, confirming that the primer set developed in this study is specific for porcine material (data not shown). However, non-specific amplifications were observed for the real-time PCR assays depending on the master mix and platform. A recent study revealed that non-specific amplifications during real-time PCR assays may occur after 27 cycles depending on the master mix and thermal cycling conditions, which is earlier than the reported threshold of 34 cycles (Ruiz-Villalba et al., 2017). In the current study, non-specific amplifications occurred with C<sub>T</sub> values ranging from 32.75 to 34.91, which are consistent with the data from an earlier study. However, the reason why certain master mixes result in non-specific amplifications remains unknown.

Table 3. Comparison of the specificity of the seven commercial master mixes for detecting porcine DNA by real-time PCR

N. C.	DI (C	C <sub>T</sub> values <sup>1,2)</sup> of animal species					
Manufacturer	Platform -	Pig	Human	Dog	Cow	Chicken	
Applied Biosystems	StepOnePlus	$17.77 \pm 0.53$	ND	ND	ND	ND	
	CFX Connect	$17.38 \pm 0.23$	ND	ND	ND	ND	
CancerROP	StepOnePlus	18.54±0.09	33.21±0.54	ND	34.28 (2/3)	ND	
	CFX Connect	$15.72\pm0.44$	34.11 (2/3)	ND	ND	ND	
Invitrogen	StepOnePlus	17.01±0.16	ND	ND	ND	ND	
	CFX Connect	$17.60\pm0.34$	ND	ND	34.28 (2/3)	ND	
Kogene Biotech	StepOnePlus	15.35±0.37	ND	ND	ND	ND	
	CFX Connect	$17.40 \pm 0.45$	ND	ND	ND	ND	
New England Biolabs	StepOnePlus	17.27±0.04	ND	ND	ND	ND	
	CFX Connect	17.61±0.11	$33.56 \pm 0.83$	ND	34.91 (2/3)	ND	
Qiagen	StepOnePlus	17.19±0.23	33.96 (2/3)	ND	ND	ND	
	CFX Connect	$17.08\pm0.26$	34.34 (2/3)	$32.75 \pm 0.73$	33.92±0.33	ND	
Takara	StepOnePlus	18.69±0.37	ND	ND	ND	ND	
	CFX Connect	$16.68 \pm 0.18$	ND	ND	ND	ND	

<sup>1)</sup> Average C<sub>T</sub> values (mean±SD) for 50 ng/µL DNA samples of five species analyzed in triplicate.

## Application of the real-time PCR assay for detecting porcine material

Following a comparison of seven master mixes and platforms based on the LOD, PCR efficiency, and specificity, we selected PowerAmp Real-time PCR Master Mix II and the CFX Connect platform as the best combination for detecting porcine DNA. Pork is one of the most widely used meats in processed meat products. The presence of pork in processed foods may cause allergic reactions in some sensitive people, and it may be used for food adulteration (Soares et al., 2010; Tanabe et al., 2007). Furthermore, it is banned in halal foods.

The real-time PCR assay developed in this study was used to determine the presence or absence of porcine material in various processed foods. The analysis of 10 commercial processed meat products showed positive reactions for all pork-containing products (dumpling, ham, pork cutlet, and sausage) and negative reactions for all beef-containing products (beef curry, beef stock, and meatballs) and chicken-containing products (chicken teriyaki, chicken sausage, and chicken stock) (Table 4). This result was consistent with the ingredients listed on the label of the commercial products. Therefore, the porcine-specific real-time PCR assay could be used for the identification and detection of hidden allergens and food adulterants in processed foods and the analysis of halal foods.

#### Comparison of costs and times

The seven master mixes included in this study were compared in terms of their costs and times required for the real-time PCR (Table 5). The cost per reaction of the seven master mixes ranged from US\$ 0.33 to 1.26. The most and least expensive master mixes were from Invitrogen and CancerROP, respectively. On the basis of the manufacturer's recommended protocols, the total real-time PCR run-times ranged from 10.75 to 57.50 min, with the longest and shortest run-times associated with the

 $<sup>^{2)}</sup>$  Numbers in parentheses indicates the number of times the  $C_T$  value was determined in three measurements, the average  $C_T$  value is presented. ND, not detected.

Table 4. Porcine-specific real-time PCR assay<sup>1)</sup> of commercial processed meat products

Processed meat products	Labeled meat ingredients	C <sub>T</sub> values <sup>2)</sup>
Dumpling	Pork	21.71±0.73
Ham	Pork	19.87±0.64
Pork cutlet	Pork	$20.06 \pm 0.42$
Sausage	Pork	18.62±0.23
Beef curry	Beef	ND
Beef stock	Beef	ND
Meatballs	Beef	ND
Chicken teriyaki	Chicken	ND
Chicken sausage	Chicken	ND
Chicken stock	Chicken	ND

<sup>1)</sup> PowerAmp Real-time PCR Master Mix II and the CFX connect platform were used for porcine-specific real-time PCR assay.

Table 5. Comparison of the time and cost associated with the seven commercial master mixes included in this study

Manufacturer	Master mix	Cost/reaction (US\$) <sup>1)</sup>	Total time (min) <sup>2)</sup>
Applied Biosystems	TaqMan Universal PCR Master Mix	1.13	38.25
CancerROP	MG 2X qPCR MasterMix (TaqMan) with ROX	0.33	57.50
Invitrogen	Express qPCR Supermix Universal	1.26	47.75
Kogene Biotech	PowerAmp Real-time PCR Master Mix II	0.59	55.75
New England Biolabs	Luna Universal Probe qPCR Master Mix	0.47	27.25
Qiagen	QuantiNova Probe PCR Kit	0.85	10.75
Takara	Premix Ex Taq (Probe qPCR), ROX plus	0.49	12.58

<sup>&</sup>lt;sup>1)</sup> Calculated based on the cost for each commercial master mix in spring 2020.

CancerROP and Qiagen master mixes, respectively.

#### **Conclusion**

Seven commercial master mixes were evaluated to optimize the real-time PCR assay conditions for detecting porcine materials. Additionally, LODs, PCR efficiency, and specificity, as well as the total times and costs of each master mix were analyzed using two real-time PCR platforms. Consequently, a sensitive (LOD of 0.5 to 5 pg/rxn) and efficient (84.96% to 108.80%) real-time PCR-based porcine detection system was developed. However, the LODs and PCR efficiencies varied depending on the master mixes and platforms. Moreover, a specificity test involving four animal species unrelated to pig revealed that non-specific amplifications were not observed on both platforms for only three master mixes from Applied Biosystems, Kogene Biotech, and Takara. These results prove that real-time PCR assays can be influenced by the master mixes and platform. The analysis of commercial processed meat products using the porcine-specific real-time PCR assay showed results that were consistent with the ingredients listed on the label. The result of this comparative study may be useful

 $<sup>^{2)}</sup>$  Average  $\dot{C}_T$  values (mean  $\pm SD$ ) for 50 ng/ $\mu L$  DNA samples of 10 processed meat products analyzed in triplicate. ND, not detected.

<sup>&</sup>lt;sup>2)</sup> Summing of the cycling times in the manufacturer's protocol.

for optimizing porcine detection system based on a real-time PCR assay.

## **Conflicts of Interest**

The authors declare no potential conflicts of interest.

# **Acknowledgements**

This work was supported by National Research Foundation of Korea (NRF-2018R1D1A1A02086040) and Dongguk University (S-2019-G0041-00036), Korea.

#### **Author Contributions**

Conceptualization: Kang SJ, Hong KW. Data curation: Kang SJ, Jang CS. Formal analysis: Kang SJ, Jang CS, Son JM. Methodology: Kang SJ, Son JM. Writing - original draft: Kang SJ, Hong KW. Writing - review & editing: Kang SJ, Jang CS, Son JM, Hong KW.

## **Ethics Approval**

This article does not require IRB/IACUC approval because there are no human and animal participants.

#### References

- Adams PS. 2006. Data analysis and reporting. In Real-time PCR. Dorak MT (ed). Taylor and Francis, New York, NY, USA. pp 39-62.
- Ali ME, Hashim U, Mustafa S, Che Man YB, Dhahi TS, Kashif M, Uddin MK, Abd Hamid SB. 2012. Analysis of pork adulteration in commercial meatballs targeting porcine-specific mitochondrial cytochrome b gene by TaqMan probe real-time polymerase chain reaction. Meat Sci 91:454-459.
- Amaral JS, Santos G, Oliveira MBPP, Mafra I. 2017. Quantitative detection of pork meat by EvaGreen real-time PCR to assess the authenticity of processed meat products. Food Control 72:53-61.
- Arslan A, Ilhak I, Calicioglu M, Karahan M. 2005. Identification of meats using random amplified polymorphic DNA (RAPD) technique. J Muscle Foods 16:37-45.
- Bustin SA. 2004. A-Z of quantitative PCR. International University Line (IUL), La Jolla, CA, USA. pp 87-112.
- Buzard GS, Baker D, Wolcott MJ, Norwood DA, Dauphin LA. 2012. Multi-platform comparison of ten commercial master mixes for probe-based real-time polymerase chain reaction detection of bioterrorism threat agents for surge preparedness. Forensic Sci Int 223:292-297.
- Eischeid AC, Kasko SM. 2015. Quantitative multiplex real-time PCR assay for shrimp allergen: Comparison of commercial master mixes and PCR platforms in rapid cycling. J Food Prot 78:230-234.
- El Sheikha AF, Mokhtar NFK, Amie C, Lamasudin DU, Isa NM, Mustafa S. 2017. Authentication technologies using DNA-based approaches for meats and halal meats determination. Food Biotechnol 31:281-315.
- Hedman J, Nordgaard A, Dufva C, Rasmusson B, Ansell R, Rådström P. 2010. Synergy between DNA polymerases increases

- polymerase chain reaction inhibitor tolerance in forensic DNA analysis. Anal Biochem 405:192-200.
- Hedman J, Nordgaard A, Rasmusson B, Ansell R, Rådström P. 2009. Improved forensic DNA analysis through the use of alternative DNA polymerases and statistical modeling of DNA profiles. Biotechniques 47:951-958.
- Hossain MA, Ali ME, Abd Hamid SB, Asing, Mustafa S, Mohd Desa MN, Zaidul ISM. 2016. Double gene targeting multiplex polymerase chain reaction-restriction fragment length polymorphism assay discriminates beef, buffalo, and pork substitution in frankfurter products. J Agric Food Chem 64:6343-6354.
- Kane DE, Hellberg RS. 2016. Identification of species in ground meat products sold on the U.S. commercial market using DNA-based methods. Food Control 59:158-163.
- Kesmen Z, Gulluce A, Sahin F, Yetim H. 2009. Identification of meat species by TaqMan-based real-time PCR assay. Meat Sci 82:444-449.
- Kim M, Yoo I, Lee SY, Hong Y, Kim HY. 2016. Quantitative detection of pork in commercial meat products by TaqMan<sup>®</sup> real-time PCR assay targeting the mitochondrial D-loop region. Food Chem 210:102-106.
- Kreader CA. 1996. Relief of amplification inhibition in PCR with bovine serum albumin or T4 gene 32 protein. Appl Environ Microbiol 62:1102-1106.
- Mane BG, Tanwar VK, Girish PS, Sharma D, Dixit VP. 2008. RAPD markers for differentiation of meat species. Indian J Vet Res 17:9-13.
- Mohamad AN, Mustafa S, Mokhtar NFK, El Sheikha AF. 2018. Molecular beacon-based real-time PCR method for detection of porcine DNA in gelatin and gelatin capsules. J Sci Food Agric 98:4570-4577.
- Naaum AM, Shehata HR, Chen S, Li J, Tabujara N, Awmack D, Lutze-Wallace C, Hanner R. 2018. Complementary molecular methods detect undeclared species in sausage products at retail markets in Canada. Food Control 84:339-344.
- Rahmati S, Julkapli NM, Yehye WA, Basirun WJ. 2016. Identification of meat origin in food products: A review. Food Control 68:379-390.
- Rodríguez MA, García T, González I, Hernández PE, Martín R. 2005. TaqMan real-time PCR for the detection and quantitation of pork in meat mixtures. Meat Sci 70:113-120.
- Ruiz-Villalba A, van Pelt-Verkuil E, Gunst QD, Ruijter JM, van den Hoff MJ. 2017. Amplification of nonspecific products in quantitative polymerase chain reactions (qPCR). Biomol Detect Quantif 14:7-18.
- Sakai Y, Kotoura S, Yano T, Kurihara T, Uchida K, Miake K, Akiyama H, Tanabe S. 2011. Quantification of pork, chicken and beef by using a novel reference molecule. Biosci Biotech Biochem 75:1639-1643.
- Soares S, Amaral JS, Mafra I, Oliveira MBP. 2010. Quantitative detection of poultry meat adulteration with pork by a duplex PCR assay. Meat Sci 85:531-536.
- Sohni Y, Kanjilal S, Kapur V. 2008. Performance evaluation of five commercial real-time PCR reagent systems using TaqMan assays for *B. anthracis* detection. Clin Biochem 41:640-644.
- Stephens KW, Hutchins RJ, Dauphin LA. 2010. Cross-platform evaluation of commercial real-time reverse transcription PCR master mix kits using a quantitative 5'nuclease assay for Ebola virus. Mol Cell Probes 24:370-375.
- Tanabe S, Miyauchi E, Muneshige A, MIO K, Sato C, Sato M. 2007. PCR method of detecting pork in foods for verifying allergen labeling and for identifying hidden pork ingredients in processed foods. Biosci Biotechnol Biochem 71:1663-1667.
- Wolffs P, Grage H, Hagberg O, Rådström P. 2004. Impact of DNA polymerases and their buffer systems on quantitative real-time PCR. J Clin Microbiol 42:408-411.
- Zia Q, Alawami M, Mokhtar NFK, Nhari RMHR, Hanish I. 2020. Current analytical methods for porcine identification in meat and meat products. Food Chem 324:126664.